


Product Specification

SLEE medical



Product:	Cellspin I
Art.No.:	15500001
Commodity Code:	84219100

Highlights

Flexible choice and handling of rotors

- Cyto centrifugation open or closed (12 loading positions)
- Easily removable rotor for outside loading
- Wide variety of rotors available (sediment, cyto centrifuge)

Competitively priced consumables

- Single and double cellfunnels (reusable or disposable)
- Large volume Ecofunnels (reusable or disposable)
- Absorption pads
- Cyto centrifugation Slides

Easy and safe handling

- Convenient preparation
- Slide clips applicable to reusable and disposable funnels
- Slide label readable at any time
- Indication of the RCF during centrifugation

Standard Accessories

- Basic instrument
- Cellcliprotor, closed (12 positions)
- 12x Cellclip (stainless steel)
- 12x Cellfunnel reusable
- Filter cards for Cellfunnel reusable; #39910001 Standard (200 St.)
- Operation manual
- Mains cable

Specifications

Power supply:	208 – 240 V 50 – 60 Hz
Power Consumption:	300 VA
Chamber volume:	Cellfunnel Single: 0.5 ml per chamber Cellfunnel Double: 0.5 ml per chamber ECOfunnel: 5 ml
Rotor Slide Capacity:	Cellcliprotor; Up to 12 microslides
Max. RPM (speed):	With cellcliprotor: 2.000 min ⁻¹ Centrifuge with applicable rotor: 6.000 min ⁻¹
Running time:	1 – 99 minutes, ∞ continuous run, short cycle mode (impulse key)
Safety:	Imbalance switch-off Automatic rotor recognition Lid locking and holding
Dimensions [W x D x H]:	366 mm x 430 mm x 257 mm
Weight [w/o accessories]:	24 kg

Optional Accessories

Cellfunnel reusable (12 pcs.), Art.No.: #39910000

Cellfunnel disposable (100 pcs.) , Art.No.: #39910004

Cellfunnel disposable (500 pcs.) , Art.No.: #39910005

Filtercards for Cellfunnel reusable (200 pcs.) , Art.No.: #39910001

Filtercards for Cellfunnel reusable, Thin, for samples less 0.4 ml (200 pcs.) , Art.No.: #39910006

Cellslides Single circle uncoated (100 pcs.) , Art.No.: #39910002

Cellslides Single circle coated (100 pcs.) , Art.No.: #39910003

Double Cellfunnel reusable (12 pcs.) , Art.No.: #39910010

Double Cellfunnel disposable, offset (100 pcs.) , Art.No.: #39910014

Double Cellfunnel disposable, offset (500 pcs.) , Art.No.: #39910015

Filtercards for Double Cellfunnel reusable (500 pcs.) , Art.No.: #39910011

Cellslides Double circle uncoated (100 pcs.) , Art.No.: #39910012

Cellslides Double circle coated (100 pcs.) , Art.No.: #39910013

Cellslides Double circle, offset (100 pcs.) , Art.No.: #39910016

ECOfunnel reusable (12 pcs.) , Art.No.: #39910020

ECOfunnel disposable (100 pcs.) , Art.No.: #39910023

ECOfunnel disposable (500 pcs.) , Art.No.: #39910024

ECOseal for ECOfunnel (100 pcs.) , Art.No.: #39910021

ECOslides coated (100 pcs.) , Art.No.: #39910022

Cellclip for Cellcliprotor (1 pc.) , Art.No.: #39910031

SLEE medical GmbH

Lise-Meitner-Strasse 11
55129 Mainz · Germany

Tel. +49 6131 95871-0
Fax +49 6131 95871-722

mail@slee.de
www.slee.de