

# Product Specification

**SLEE** medical



Product:	slidetec <sub>HEAT</sub>
Art.No.:	13000510 (230 V) 13000511 (115 V)
Commodity Code:	84191900

## Highlights

- Histology slide warmer
- Staged slide warming area for 30 glass microslides
- Integrated overheat security and safety thermostat
- Digital temperature display
- Temperature selection up to 60 °C
- Programmable working hours

## Standard Accessories

- Basic instrument
- Histology slide warmer (staged, aluminum)
- Operation manual
- Mains cable

# SPECIFICATIONS

Temperature range	up to 60 °C for slide warmer
Slide warmer capacity	30 microslides
Control	microprocessor, programming of temperature, working day, working time, time and date
Dimensions [W x D x H]	350 x 370 x 230 mm
Weight (w/o accessories)	15 kg
Electrical parameters	Art.No.: 13000510 230 V; 50/60 Hz, Pmax: 300 W Art.No.: 13000511 115 V; 50/60 Hz, Pmax: 300 W

SLEE medical GmbH

Lise-Meitner-Strasse 11  
55129 Mainz · Germany

Tel. +49 6131 95871-0  
Fax +49 6131 95871-722

[mail@slee.de](mailto:mail@slee.de)  
[www.slee.de](http://www.slee.de)