

# Product Specification

**SLEE** medical



Product:	MPS/P1
Art.No.:	10190000 (230 V) 10192000 (115 V)
Commodity Code:	90189084

## Highlights

- Manual lever and foot switch for paraffin dispenser
- Paraffin container 3.3 litres
- Cooling spot
- Two heating trays for 80 cassettes and 200 embedding moulds.
- Ergonomic working height
- Working area with LED light

## Standard Accessories

- Paraffin Dispensing Module MPS/P1
- Magnifying glass
- Paraffin collection tray
- Working area illumination (Light)
- Foot switch for paraffin dispenser
- Integrated cassette bath and mold warmer
- Operation manual
- Mains cable

## SPECIFICATIONS

Nominal power:	230 V AC +-10% (Art.No. 10190000) 115 V AC +-10% (Art.No. 10192000)
Nominal frequency:	50 / 60 Hz
Power:	650 W
Protective class:	(1) I
Fuses:	2x T 4 A (Art.No. 10190000) 2x T 10 A (Art.No. 10192000)
Pollution degree:	(1) 2
Overvoltage protection category:	II
Max. Heat emission:	650 J/s
Relative humidity:	max. 80%, non condensing
Humidity during storage:	< 8%
Dimensions [Wx D x H]:	450 mm x 575 mm x 295 mm
Weight [w/o accessories]:	23 kg
Paraffin container:	3.3 litres
Heated cassette/mould container:	approx. 80 cassettes approx. 200 molds 30 – 80 °C programmable
Forceps holders:	2 x 3
Illumination:	LED
Work area:	450 * 300 mm
Cooling spot:	Peltier, 50 * 50 mm

Hints for installation < 2000 m NN, 15 cm distance from backside to wall.

## Optional Accessories

Heated Forceps (2.5 mm), Art.No.: #10196000

Heated Forceps (1.25 mm), Art.No.: #10196001

Anatomic Forceps, 130 mm, stainless steel (pointed), Art.No.: #30030090

Anatomic Forceps, 130 mm, stainless steel (blunt), Art.No.: #30030091

Base moulds (stainless steel) 7 mm x 7 mm x 5 mm, Art.No.: #30030050

Base moulds (stainless steel) 30 mm x 24 mm x 5 mm, Art.No.: #30030051

Base moulds (stainless steel) 15 mm x 15 mm x 5 mm, Art.No.: #30030052

Base moulds (stainless steel) 24 mm x 24 mm x 5 mm, Art.No.: #30030053

Base moulds (stainless steel) 37 mm x 24 mm x 5 mm, Art.No.: #30030054

Base moulds (stainless steel) 33 mm x 24 mm x 12 mm, Art.No.: #30030056

Metal lids for embedding cassettes, Art.No.: #30050100

Paraffin (bag, 20 kg, melting point: 56 – 58 °C), Art.No.: #30010002

SLEE medical GmbH

Lise-Meitner-Strasse 11  
55129 Mainz · Germany

Tel. +49 6131 95871-0  
Fax +49 6131 95871-722

[mail@slee.de](mailto:mail@slee.de)  
[www.slee.de](http://www.slee.de)