

# Product Specification

**SLEE** medical



Product:	MPS/P2
Art.No.:	10185000
Commodity Code:	90189084

## Highlights

- Spacious cooling plate
- Adjustable magnifier
- Ergonomic hand rest
- Electrically heated forceps
- Programmable working hours
- Easily adjustable forceps holders
- 5 liters paraffin container
- Spacious and heated paraffin collection tray
- Working area with adjustable LED-illumination and integrated cold spot

## Standard Accessories

- Paraffin Dispensing Module MPS/P
- Magnifying glass
- Paraffin collection tray
- Working area illumination (Light)
- Foot switch for paraffin dispenser
- Prewarming module MPS/W or Prewarming Module MPS/WH
- Cooling module MPS/C
- Operation manual
- 3x Mains cable

	MPS/P		MPS/W		MPS/AWH		MPS/C		MPS/CX	
	115 V	230 V	115 V	230 V	115 V	230 V	115 V	230 V	115 V	230 V
Nominal supply voltage	115 V AC +/- 10%	230 V AC +/- 10%	115 V AC +/- 10%	230 V AC +/- 10%	115 V AC +/- 10%	230 V AC +/- 10%	115 V AC +/- 10%	230 V AC +/- 10%	115 V AC +/- 10%	230 V AC +/- 10%
Nominal frequency	60 Hz	50 / 60 Hz	60 Hz	50 / 60 Hz	60 Hz	50 / 60 Hz	60 Hz	50 / 60 Hz	60 Hz	50 / 60 Hz
Power draw	340 VA		220 VA		220 VA		110 VA		150 VA	
Protective Class (1)	I									
Power fuses	2*T2,5A	2*T2,5A	2*T2,5A	2*T2,5A	2*T2,5A	2*T2,5A	2*T6,3A	2*T4 A	2*T6,3A	2*T4A
Pollution degree (1)	2									
Overvoltage installation category	II									
Maximum heat emission	340 J/s		220 J/s		220J/s		110 J/s		150 J/s	
Operating temperature range	+ 10 bis + 35 °C									
Temperature range during storage	+ 5 bis + 55 °C									
Relative humidity	Max. 80 %, non condensing									
Humidity during storage	< 8%									
Dimensions (L*W*H) / mm	450* 575*295		300* 575* 295		300* 575* 295		400* 575* 295		600* 575* 295	
weight / kg	23		13		13		24		30	
Paraffin tank	5 l									
Forceps holders	2 * 4									
Light	LED									
Temperature range cool plate							-5 bis -15°C programmable		-5 bis -15°C programmable	
Cooling gas							R 134 a		R 134 a	
Work area							200 * 300 mm 80 cassettes, 250 moulds		200 * 450 mm 120 cassettes, 375 moulds	
Cassette bath, storage for moulds			Approx. 150 cassettes, approx. 500 moulds 30 – 80°C		Approx. 250 cassettes, approx. 750 moulds 30 – 80°C					
Installation requirements	< 2.000 m NN, 15 cm distance to wall									

(1) According to IEC 1010, EN 61010

## Optional Accessories

Heated Forceps (2.5 mm), Art.No.: #10196000

Heated Forceps (1.25 mm), Art.No.: #10196001

Anatomic Forceps, 130 mm, stainless steel (pointed), Art.No.: #30030090

Anatomic Forceps, 130 mm, stainless steel (blunt), Art.No.: #30030091

Base moulds (stainless steel) 7 mm x 7 mm x 5 mm, Art.No.: #30030050

Base moulds (stainless steel) 30 mm x 24 mm x 5 mm, Art.No.: #30030051

Base moulds (stainless steel) 15 mm x 15 mm x 5 mm, Art.No.: #30030052

Base moulds (stainless steel) 24 mm x 24 mm x 5 mm, Art.No.: #30030053

Base moulds (stainless steel) 37 mm x 24 mm x 5 mm, Art.No.: #30030054

Base moulds (stainless steel) 33 mm x 24 mm x 12 mm, Art.No.: #30030056

Metal lids for embedding cassettes, Art.No.: #30050100

Paraffin (bag, 25 kg, melting point: 56 – 58 °C), Art.No.: #30010002

SLEE medical GmbH

Lise-Meitner-Strasse 11  
55129 Mainz · Germany

Tel. +49 6131 95871-0  
Fax +49 6131 95871-722

[mail@slee.de](mailto:mail@slee.de)  
[www.slee.de](http://www.slee.de)