

## Product Specification

**SLEE** medical



Product:	slidetec <sub>WATER/HEAT</sub>
Art.No.:	13000500 (230 V) 13000501 (115 V)
Commodity Code:	84191900

### Highlights

- Histology water bath and slide warmer combination
- Separable glass water bath
- Integrated water bath illumination (LED, white/blue)
- Integrated overheat security and safety thermostat
- Digital temperature display
- Temperature selection up to 60 °C
- Programmable working hours

### Standard Accessories

- Basic instrument
- Histology water bath (glass, LED illumination)
- Slide warmer (stage)
- Operation manual
- Mains cable

# SPECIFICATIONS

Temperature range	up to 60 °C for water bath up to 60 °C for slide warmer
Water bath capacity	2.0 litres
Water bath illumination	LED (white/blue)
Slide warmer capacity	30 slides
Control	microprocessor, programming of temperature, working day, working time, time and date
Dimensions [W x D x H]	560 x 370 x 230 mm
Weight (w/o accessories)	18 kg
Electrical parameters	Art.No.: 13000500 230 V; 50/60 Hz, Pmax: 600 W Art.No.: 13000501 115 V; 50/60 Hz, Pmax: 600 W

SLEE medical GmbH

Lise-Meitner-Strasse 11  
55129 Mainz · Germany

Tel. +49 6131 95871-0  
Fax +49 6131 95871-722

mail@slee.de  
www.slee.de